

# NOBAMED Paul Danz AG

# **RUDASTRIP®**

### **REF 066006**

# **Product Description, Intended use, Application**

Reliably adhesive white-coloured wound suture strips, size 6 mm x 100 mm, packed per 10 pieces in sterile packaging. The wound suture strips are used for an atraumatic primary closure of small, smooth, superficial lacerations or cuts, surgical incisions, as a complement to/relief for subcutaneous and intra-cutaneous sutures or for wounds that have been closed with staples, for the fixation of skin grafts, and if staples and sutures are prematurely removed. RUDASTRIP® are for single use only.

### Composition

Polyamide, polyester, polyacrylate adhesive

### **Contra-Indications**

Wounds under tension, infected wounds, deeper wounds, and heavily bleeding wounds should not be closed with wound suture strips.

The product should not be used in the case of a known allergy against the material.

### **Note**

RUDASTRIP® is not intended for the fixation of e.g., dialysis cannulae.

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

# Normative and Regulative Requirements, Common Standards

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

Sterilization of the product complies with DIN EN ISO 11135.

The product does not contain dangerous toxic substances according to REACH.

## **Packaging**

Primary packaging: pa

paper-filmpackaging

Secondary packaging: folding box

made of cellulose

Tertiary packaging: carton made

of cellulose

### Symbols used in labelling

Explanation at www.nobamed.com











**€** 

Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417

### Storage and Transport

Dry and dustfree

### Sterile device

Before using a sterile product, visually inspect the packaging to ensure that it is intact.

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.