

NOBAMED Paul Danz AG

RUDASTRIP®

REF 066012

Product Description, Intended use, Application

Reliably adhesive white-coloured wound suture strips, size 13 mm x 100 mm, packed per 6 pieces in sterile packaging. The wound suture strips are used for an atraumatic primary closure of small, smooth, superficial lacerations or cuts, surgical incisions, as a complement to/relief for subcutaneous and intra-cutaneous sutures or for wounds that have been closed with staples, for the fixation of skin grafts, and if staples and sutures are prematurely removed. RUDASTRIP® are for single use only.

Composition

Polyamide, polyester, polyacrylate adhesive

Contra-Indications

Wounds under tension, infected wounds, deeper wounds, and heavily bleeding wounds should not be closed with wound suture strips.

The product should not be used in the case of a known allergy against the material.

Note

RUDASTRIP® is not intended for the fixation of e.g., dialysis cannulae.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulative Requirements, Common Standards

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

Sterilization of the product complies with DIN EN ISO 11135.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging: pap

paper-filmpackaging

Secondary packaging: folding box

made of

cellulose

Tertiary packaging:

carton made of cellulose

Symbols used in labelling

Explanation at www.nobamed.com











€

Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417

Storage and Transport

Dry and dustfree

Sterile device

Before using a sterile product, visually inspect the packaging to ensure that it is intact.

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.