



NOBASET®

Dialysekathetertaschen- verband

REF 075125 (ab LOT 35887)

Content of this procedure pack:

1 x protection pouch for dialysis catheters,

NOBAMED Paul Danz AG, 

1 x Rudavlies 12,5 x 2,5 cm, NOBAMED

Paul Danz AG, 

1 x Rudaven-Inject 6 x 9 cm, NOBAMED

Paul Danz AG, 

Product Description, Intended use, Application

The dialysis catheter pouch dressing is used to fix dialysis catheters quickly and safely and to protect them against mechanical stress and contamination. By RUDAVEN®-inject the catheter is fixed at the exit site. The catheter is placed in the pouch which is then fixed to the skin using RUDAVLIES®.

Composition

Protection pouch: Polyester, elastane, RUDAVLIES®: Polyester, polyacrylate adhesive, RUDAVEN®-Inject: Viscose, polyester, polyethylene, polyacrylate adhesive

Contraindication:

Do not use on infected wounds. Not suitable for wounds that bleed or exude heavily. The product should not be used in the case of a known allergy against the material.



This product is for single and temporary use on one patient and must be applied by specially trained medical personnel. Any information for use that is available for individual products must be observed!

Normative and Legal Requirements

Procedure Pack according to Art. 22 of the Medical Device Regulation (EU) 2017/745, risk class I.

The compatibility of the customized medical devices or other products has been assessed with due regard for the indications of the respective manufacturers, and operations have been performed in accordance with the aforementioned indications.

The manufacturers' intended use is not changed.

The systems and procedure packs are packed in compliance with a validated procedure.

All activity is properly monitored and controlled as part of our QM system, in accordance with DIN EN ISO 13485:2016.

A validated sterilization process is used for sterile procedure packs, taking into account the manufacturers' instructions, so that product functionality and product safety of the components are maintained after a system or procedure pack has been compiled and sterilized.

Sterilization and packaging are subject to supervision of our Notified Body TÜV Süd PS 0123 pursuant to Annex XI, Part A.

The product does not contain dangerous toxic substances according to REACH. It has CE marking and DIN EN ISO 15223-1 labels on all its packaging

Storage and transport

In a dry and dust-free environment, protect from sunlight.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.

The product bears following labels and symbols Explanation at www.nobamed.com

