



## RUDALASTIK®

REF 095506

### Product Description, Intended use, Application

The longitudinal elastic plaster bandage (width: 6 cm, unstretched length: 2.5 m) is coated with a hypoallergenic adhesive. The adhesive layer does not impair the open-pored structure of the bandage. RUDALASTIK® stabilizes and compresses muscles, tendons and joints and is in most cases used together with RUDATAPE® for the application of functional bandages. In phlebology, RUDALASTIK® is used as a compression in cases of venous insufficiency.

### Composition

Cotton, polyacrylate adhesive

### Contraindications

The product should not be used in the case of a known allergy against the material.

### Residual Risks, Undesirable Side Effects, Warnings, Notes

In the case of a medical compression therapy the contra-indications and risks mentioned in the relevant medical guidelines must be observed.

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Regulatory Requirements, Common Standards

Medical Device according to *MDD 93/42/EEC*, *MDR (EU) 2017/745*.

The product does not contain dangerous toxic substances according to REACH.

### Packaging

Primary packaging: folding box made of cellulose

Secondary packaging: carton made of cellulose

### Symbols used in labelling

Explanation at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417

### Storage and Transport

Dry and dustfree, protect from sunlight.

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.