



NOBAMED Paul Danz AG

NOBADERM®

REF 160606

Product Description and Purpose

The sterile, transparent film dressing, 6 cm x 7 cm net dressing area, consists of polyurethane and is coated with a hypoallergenic polyacrylate adhesive. NOBADERM® is suitable for moist wound treatment. Being applied as a primary dressing it protects weakly secreting wounds from the penetration of germs. It is applied on wounds in the epithelisation phase and on postoperative wounds. NOBADERM® is also applied as a secondary dressing being excellently suitable for moist wound treatment in combination with for example alginates (NOBAALGIN®) and hydrogels (NOBAGEL®). The semipermeable product is waterproof, but water vapour permeable (MVTR ~ 900-1000 g/sqm/24h). NOBADERM® is individually packed.

Contra-Indications

May not be used on infected wounds. The dressing is not suitable as a primary dressing for heavily bleeding or moderate exsuding wounds. May not be used in case of allergies against one of the ingredients listed below.

Composition

Polyurethane film, polyacrylate adhesive

Normative and Legal Requirements

NOBADERM® is a medical product according to the MPG, the directive 93/42/EEC and the regulation MDR (EU 2017/745) and is classified as class IIa, rule 4.

Sterilization of the product complies with DIN EN 11135.

The product does not contain dangerous toxic substances according to REACH.

It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

Packaging

Primary packaging:	paper- film packaging
Secondary packaging:	folding box made of cellulose
Tertiary packaging:	carton made of cellulose

Storage

To be stored in a dry and dust-free environment at room temperature, protected from heat and sunlight

The product bears following labels and symbols:

