

**NOBAMED** Paul Danz AG

# **NOBADERM®**

# **REF 161012**

# Product Description, Intended use, Application

NOBADERM<sup>®</sup> (10 cm x 12 cm) is a sterile, transparent film dressing. It is made of polyurethane and is equipped with an adhesive coating based on polyacrylate. NOBADERM<sup>®</sup> as primary wound dressing protects wounds from the penetration of germs and serves as catheter fixation.

# Composition

Polyurethane film, adhesive based on polyacrylate

# **Contraindications**

Do not use NOBADERM<sup>®</sup> for infected wounds. Not suited as primary wound dressing for heavily bleeding or highly exuding wounds. Do not use if there is a known allergy to one of the ingredients of NOBADERM<sup>®</sup>.

### **Incident reporting**

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## Normative and Regulative Requirements, Common Standards

Medical Device according to MDD 93/42/EEC and MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

Sterilization of the product complies with DIN EN ISO 11135.

# Packaging

Primary:	paper-film packaging
Secondary:	folding box made of cellulose
Tertiary:	carton made of cellulose

# Symbols used in labelling

Explanations at <u>www.nobamed.com</u>



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417.

# **Storage and Transport**

Expiry date on all packaging levels, protect against moisture, sunlight and heat, store in a dust-free environment at room temperature. The product must no longer be used after the expiry date has passed.

### **Sterile device**

Contents sterile only if package is intact. If the primary package is damaged, the product may no longer be used.

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.