



**NOBAMED Paul Danz AG**

## Woundhook cover

**REF 324022** (as of LOT 137922)

### Product Description, Intended use, Application

Disposable single-layer white covers sewn of medically bleached knitted cotton fabric, 3 cm x 22 cm, packed per 20 pcs, with worked in X-ray thread and with bands. The bands are used to fix the covers at the wound hook. The covers are put on the wound hooks before sterilization, to protect the instruments. The covers may be left on the wound hooks during their use and will thus reduce light reflections caused by the metal due to the operation lamp. At the same time sensitive tissue will be protected against the edges of the wound hooks. Before use on open wounds, the product must be sterilized in accordance with a validated sterilization method (Moist Heat in accordance with DIN EN ISO 17665, Radiation in accordance with DIN EN ISO 11137, or Ethylene Oxide in accordance with DIN EN ISO 11135).

### Composition

cotton, x-ray thread made of barium sulfate, PVC, phthalate-free plasticizer

### Note

Sterilize before use

### Contraindications

Allergy caused by the material

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Regulatory Requirements, Common Standards

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

### Packaging

Primary packaging: PE bag

Secondary packaging: carton made of cellulose

### Symbols used in labelling

Explanation at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417

### Storage and Transport

Dry and dustfree, protect from sunlight.

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.