

## **NOBAMED Paul Danz AG**

## NOBALINE®-Tuch

### **REF 771304**

# **Product Description, Intended use, Application**

NOBALINE®-Tuch is a 35 cm x 40 cm nonwoven dressing drape. The non-woven fabric is coated with aluminium on one side and is therefore non-adhesive and has a hydrophobic reverse side. Due to the porous structure of the aluminium layer, wound secretions can be quickly transferred to the absorbent core of the dressing cloth. In addition, the drapes are provided with elastic bandages at the corners for fixation to the patient. The dressing drape is used for sterile covering of large, sensitive wounds and injuries, as a wound dressing for open wound treatment and as an underlay to protect against germs during first aid. NOBALINE®-Tuch is intended for single use and is individually sealed.

### Composition

Polyethylene, aluminum, polyester, rayon

## **Contraindications**

The product should not be used in the case of a known allergy against the material.

### **Incident reporting**

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## Normative and Regulative Requirements, Common Standards Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

Sterilization of the product complies with DIN EN 11135.

The product does not contain dangerous toxic substances according to REACH.

## **Packaging**

Primary packaging: paper-film

packaging

Secondary packaging: PE-bag

Tertiary packaging: carton made of

cellulose

## Symbols used in labelling

Explanation at www.nobamed.com









Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417

## **Storage and Transport**

Dry and dustfree

## Sterile device

Before using a sterile product, visually inspect the packaging to ensure that it is intact.

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

## Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.

Product Data Sheet: September 13, 2022 [REV 5] replaces: November 09, 2017 [REV 4]