



**NOBAMED Paul Danz AG**

## NOBACOLLOID®

REF 790115 (as of LOT 137289)

### Product Description, Intended use, Application

NOBACOLLOID® (15 cm x 15 cm) is a hydrocolloid wound dressing.

The dressing consists of an absorbent hydrocolloid matrix and a breathable polyurethane film backing. A layer of wet gel forms once the hydrocolloid dressing gets in contact with the wound exudate, providing favourable conditions for the healing process. NOBACOLLOID® is indicated for the local treatment of leg ulcers and pressure ulcers.

### Composition

Carboxymethyl cellulose matrix, pressure sensitive adhesive, polyurethane film

### Contraindication

Do not use for infected wounds.  
Do not use the product on any patients who may be allergic to any of its ingredients.

### Precaution

The use of NOBACOLLOID® for leg ulcers does not replace the need for compression treatment (stockings or bandages) where prescribed.

The use of NOBACOLLOID® on pressure ulcers does not replace the need for normal nursing care.

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Regulative Requirements, Common Standards

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

Sterilization of the product complies with DIN EN 11137.


The product does not contain dangerous toxic substances according to REACH.

### Packaging

Primary:	foil-paper packaging
Secondary:	folding box made of cellulose
Tertiary:	carton made of cellulose

### Symbols used in labelling

Explanation at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417

### Storage and Transport

Store NOBACOLLOID® flat.  
Store with primary packing sealed.  
Keep dry, keep away from sunlight.  
NOBACOLLOID® should be stored in dry conditions at 0~25°C

**Sterile device**

Before using a sterile product, visually inspect the packaging to ensure that it is intact.

**Single use device**

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

**Disposal**

In accordance with the requirements of the applicable local legislation and with the standards for the prevention of infections.