



**NOBAMED Paul Danz AG**

# **NOBA<sup>®</sup> EXTENSIONSVERBAND (SKIN TRACTION KIT)**

without adhesive

**REF 944901**

## **Product Description, Intended Use, Application**

NOBA<sup>®</sup>-EXTENSIONSVERBAND (SKIN TRACTION KIT) **without adhesive** is a ready-to-use extension bandage for **adults**, consisting of a stirrup lined with foam, two extension strappings, an affixed string and a bandage. Both extension strappings are provided with foam material so that they do not slip out of position on the skin. The skin traction kit is used for the treatment of fractures of the thigh or lower leg bones. It is fixed to the lower leg with the enclosed bandage. Weights can be attached to the affixed strings via a suitable roller system. This allows for a continuous longitudinal traction that is exerted on the injured extremity. The skin traction kit is adequately padded in the ankle area. NOBA<sup>®</sup>-EXTENSIONSVERBAND is for single use and is individually packed.

## **Composition**

Polyester, cotton, synthetic latex, polyurethane, polypropylene, polyamide

## **Contraindications**

In the case of an enormous dislocation of the bone fragments which can only be reduced by means of a surgical intervention, and in cases of comminuted fractures.

The product should not be used in the case of a known allergy against the material.

## **Incident reporting**

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## **Normative and Regulatory Requirements, Common Standards**

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

## **Packaging**

Primary packaging: bandage in OPP-film and folding box made of cellulose

Secondary packaging: carton made of cellulose

## **Symbols used in labelling**

Explanations at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417.

## **Storage and Transport**

To be stored in a dry and dust-free environment, protected against direct sunlight

## **Single use device**

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

## **Disposal**

According to locally applicable legal regulations and standards of infection prophylaxis