



NOBAMED Paul Danz AG

Processing instructions of NOBAMED Paul Danz AG for the medical device **safety pin with steel pinhead**

Depending on the application, the product must be sterilized before use. We confirm the one-time sterilizability with the procedure mentioned under the section 'sterilization'.

Instructions	
Initial treatment at the place of use (storage and transport)	The storage and transport conditions of the non-sterile product are indicated on the package
If sterilisation is required, preparation for sterilisation	Remove and discard the packaging of the product.
Cleaning/ Disinfection	Not applicable
Drying	Not applicable
Control and Verification	The validation of the sterilization process and the product fall under the responsibility of the processor
Packaging	Suitable validated sterile barrier system
Sterilization	Validated sterilization process for example <ul style="list-style-type: none">with moist heat according to DIN EN ISO 17665 or as per requirements of the European Pharmacopoeia with 15' to 20', 121°C to 134°C, 2 bar
Storage and Transport	For the product the storage and transport conditions indicated on the original packaging are applicable.
Additional Information	Safety pins with steel pinhead are disposable products and cannot be reprocessed for product-related reasons. This does not refer to the process of first-time sterilization, if required before use.

The above-mentioned instructions have been compiled in accordance with EN ISO 17664 and have been validated as SUITABLE for the processing of Safety pins with steel pinhead. It falls under the processor's responsibility to make sure that the actual process (with the equipment, materials, and staff in the processing department) is adequate to achieve the desired results. This necessitates a validation and routine monitoring of the process on site.

Hersteller: **NOBAMED Paul Danz AG**, Höltkenstr. 1-5, 58300 Wetter (Ruhr),
Phone. (+49)02335/7609-0; Fax: (+49)02335/7609-20;
E-Mail: info@nobamed-ag.com, website: www.nobamed.com