



**NOBAMED Paul Danz AG**

## RIBOSORB®-plus

REF 001139

### Product Description, Intended use, Application

The pulp diapers (37 cm x 11 cm) are used to absorb secretion or urine e.g. in case of incontinence. They are suitable for adults and children alike. They may also serve as an inlay for absorbing exudate and transudate such as bleedings or discharge in the gynecological field. The additional use of superabsorbents in the absorbent body of the non-woven diapers increases the absorptive capacity of RIBOSORB®-plus. The outer layer consists of polypropylene non-woven.

### Composition

Cellulose, polypropylene

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Regulative Requirements, Common Standards

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

### Packaging

Primary packaging: PE bag  
Secondary packaging: Carton made of cellulose

### Symbols used in labelling

Explanations at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN 1041.

### Storage and Transport

Dry and dustfree.

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.