



## RUDAVLIES®

REF 070102

### Product Description, Intended use, Application

RUDAVLIES® is a transversely stretchable soft and white nonwoven fabric coated with hypoallergenic polyacrylate adhesive. 10 m are on a roll (width 2.5 cm). RUDAVLIES® can be cut to size individually and is used to fix primary wound dressings. In addition, RUDAVLIES® can be used for the fixation of e.g. catheters, tubes and measuring elements.

### Composition

Polyester, polyacrylate adhesive

### Contraindications

The product should not be used in the case of a known allergy against the material.

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Regulative Requirements, Common Standards

Medical Device according to *MDR (EU) 2017/745*.

The product does not contain dangerous toxic substances according to REACH.

### Packaging

Primary: folding box  
made of cellulose

Secondary: carton  
made of cellulose

### Symbols used in labelling

Explanation at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417

### Storage and Transport

Dry and dustfree

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis