



**NOBAMED Paul Danz AG**

## RUDAVLIES®

**REF 070110**

### **Product Description, Intended use, Application**

RUDAVLIES® is a transversely stretchable soft and white nonwoven fabric coated with hypoallergenic polyacrylate adhesive. 10 m are on a roll (width 10 cm). RUDAVLIES® can be cut to size individually and is used to fix primary wound dressings. In addition, RUDAVLIES® can be used for the fixation of e.g. catheters, tubes and measuring elements.

### **Composition**

Polyester, polyacrylate adhesive

### **Contraindications**

The product should not be used in the case of a known allergy against the material.

### **Incident reporting**

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### **Normative and Regulatory Requirements, Common Standards**

Medical Device according to *MDR (EU) 2017/745*.

The product does not contain dangerous toxic substances according to REACH.

### **Packaging**

Primary: folding box  
made of cellulose

Secondary: carton  
made of cellulose

### **Symbols used in labelling**

Explanation at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417

### **Storage and Transport**

Dry and dustfree

### **Single use device**

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

### **Disposal**

According to locally applicable legal regulations and standards of infection prophylaxis