

NOBAMED Paul Danz AG

# **RUDAVLIES<sup>®</sup>**

## **REF 070130**

# Product Description, Intended use, Application

RUDAVLIES<sup>®</sup> is a transversely stretchable soft and white nonwoven fabric coated with hypoallergenic polyacrylate adhesive. 10 m are on a roll (width 30 cm). RUDAVLIES<sup>®</sup> can be cut to size individually and is used to fix primary wound dressings. In addition, RUDAVLIES<sup>®</sup> can be used for the fixation of e.g. catheters, tubes and measuring elements.

### Composition

Polyester, polyacrylate adhesive

#### **Contraindications**

The product should not be used in the case of a known allergy against the material.

#### **Incident reporting**

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

#### Normative and Regulative Requirements, Common Standards Medical Device according to *MDR* (*EU*) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

# Packaging

Primary:

folding box made of cellulose

Secondary:

carton made of cellulose

Symbols used in labelling Explanation at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417

### **Storage and Transport**

Dry and dustfree

#### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

#### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis