

NOBAMED Paul Danz AG

RUDAKITAPE[®]

REF 071005

Product Description, Intended use, Application

RUDAKITAPE® (nature, 5,00 cm x 5 m) is a longitudinally elastic tape bandage that does not restrict the movement of joints. It is applied as a kinesiotape. Kinesiotaping is used to reduce pain and support muscle-ligament-joint function.

Composition

Cotton, elasthan, acrylic adhesive

Contraindications

Extensive hematomas, extensive skin damage and known allergies to any component of the product.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulative

Requirements, Common Standards Medical Device according to *MDR (EU)* 2017/745.

The product does not contain dangerous toxic substances according to REACH.

Packag	ing
Primary	:

folding box made of cellulose

Secondary:

carton made of cellulose

Symbols used in labelling Explanation at <u>www.nobamed.com</u>



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1and DIN EN ISO 20417

Storage and Transport Dry and dust free

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.