



## RUDAKITAPE®

REF 071015

### Product Description, Intended use, Application

RUDAKITAPE® (blue, 5,00 cm x 5 m) is a longitudinally elastic tape bandage that does not restrict the movement of joints. It is applied as a kinesiotape. Kinesiotaping is used to reduce pain and support muscle-ligament-joint function.

### Composition

Cotton, elasthan, acrylic adhesive

### Contraindications

Extensive hematomas, extensive skin damage and known allergies to any component of the product.

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Regulative Requirements, Common Standards

Medical Device according to *MDR (EU) 2017/745*.

The product does not contain dangerous toxic substances according to REACH.

### Packaging

Primary: folding box made of cellulose

Secondary: carton made of cellulose

### Symbols used in labelling

Explanation at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417

### Storage and Transport

Dry and dust free

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.