

# NOBAMED Paul Danz AG

## **RUDAKITAPE®**

#### **REF 071065**

# **Product Description, Intended use, Application**

RUDAKITAPE® (pink, 5,00 cm x 5 m) is a longitudinally elastic tape bandage that does not restrict the movement of joints. It is applied as a kinesiotape. Kinesiotaping is used to reduce pain and support muscle-ligament-joint function.

#### Composition

Cotton, elasthan, acrylic adhesive

#### **Contraindications**

Extensive hematomas, extensive skin damage and known allergies to any component of the product.

### **Incident reporting**

2017/745.

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## Normative and Regulative Requirements, Common Standards Medical Device according to MDR (EU)

The product does not contain dangerous toxic substances according to REACH.

## **Packaging**

Primary: folding box made of

cellulose

Secondary: carton made of cel-

lulose

#### Symbols used in labelling

Explanation at <a href="https://www.nobamed.com">www.nobamed.com</a>









Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417

## **Storage and Transport**

Dry and dust free

#### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

## Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.