

NOBAMED Paul Danz AG

RUDAVLIES®-steril

REF 072593

Product Description, Intended use, Application

The carrier material of the sterile wound dressing (size: 30 cm x 10 cm) consists of a white, elastic non-woven and is coated with a hypoallergenic polyacrylate adhesive. The wound pad does not adhere to the wound and is completely protected by the carrier material. The dressings is individually packed in sterile packaging.

Composition

Polyester, viscose, polyethylene, polyacrylate adhesive

Contraindications

The product should not be used in the case of a known allergy against the material.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulative Requirements, Common Standards

Medical Device according to MDD 93/42/EEC and MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

Sterilization of the product complies with DIN EN ISO 11135.

Packaging

Primary: paper packaging

Secondary: folding box made of

cellulose

Tertiary: carton made of cellulose

Symbols used in labelling

Explanation at www.nobamed.com

















Kennzeichnung auf allen Verpackungsstufen mit CE und nach DIN EN ISO 15223-1 und DIN EN ISO 20417.

Storage and Transport

Dry and dustfree, protect from sunlight.

Sterile device

Before using a sterile product, visually inspect the packaging to ensure that it is intact.

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.