



## RUDAVLIES®-steril

REF 072595

### Product Description, Intended use, Application

The carrier material of the sterile wound dressing (size: 35 cm x 10 cm) consists of a white, elastic non-woven and is coated with a hypoallergenic polyacrylate adhesive. The wound pad does not adhere to the wound and is completely protected by the carrier material. The dressings is individually packed in sterile packaging.

### Composition

Polyester, viscose, polyethylene, polyacrylate adhesive

### Contraindications

The product should not be used in the case of a known allergy against the material.

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Regulative Requirements, Common Standards

Medical Device according to MDD 93/42/EEC and MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

Sterilization of the product complies with DIN EN ISO 11135.

### Packaging


Primary: paper packaging

Secondary: folding box made of cellulose

Tertiary: carton made of cellulose

### Symbols used in labelling

Explanation at [www.nobamed.com](http://www.nobamed.com)



Kennzeichnung auf allen Verpackungsebenen mit CE und nach DIN EN ISO 15223-1 und DIN EN ISO 20417.

### Storage and Transport

Dry and dustfree, protect from sunlight.

### Sterile device

Before using a sterile product, visually inspect the packaging to ensure that it is intact.

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.