



## RUDAMED®-elastisch

REF 080506

### Product Description, Intended use, Application

The beige, transversely elastic adhesive dressing (6 cm x 5 m) is made of textile fabric which is provided with a well-adhering zinc oxide adhesive which has a high ageing resistance. RUDAMED®-elastisch is used to cover smaller primary healing wounds. The non-adherent wound pad absorbs blood and exudate. RUDAMED®-elastisch is packed in a convenient dispenser box with opening.

### Composition

Polyester, viscose, polyethylene, zinc oxide adhesive

### Contra indications

The product should not be used in the case of a known allergy against the material.

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Legal Requirements

Medical product according to the directive MDD 93/42/EEC and the regulation MDR (EU 2017/745).

The product does not contain dangerous toxic substances according to REACH.

### Packaging

Primary:	folding box made of cellulose
Secondary:	carton made of cellulose

### Symbols used in labelling

Explanations at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and ISO 20417.

### Storage and Transport

Dry and dustfree

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.