



## RUDAMED<sup>®</sup>-soft

REF 081504

### Product Description, Intended use, Application

The beige adhesive dressing (4 cm x 5 m) is made of an elastic, soft non-woven fabric and is provided with a hypoallergenic polyacrylate adhesive. The dressing is therefore suitable for sensitive skin. RUDAMED<sup>®</sup>-soft is used to cover smaller primary healing wounds. The non-adherent wound pad absorbs blood and exudate. RUDAMED<sup>®</sup>-soft is packed in a convenient dispenser box with opening.

### Composition

Polyester, polyethylene, viscose and polyacrylate adhesive

### Contra indications

The product should not be used in the case of a known allergy against the material.

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Legal Requirements

Medical product according to the directive MDD 93/42/EEC and the regulation MDR (EU 2017/745).

The product does not contain dangerous toxic substances according to REACH.

### Packaging

Primary:	folding box made of cellulose
Secondary:	carton made of cellulose

### Symbols used in labelling

Explanations at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and ISO 20417.

### Storage and Transport

Dry and dustfree

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.