



RUDAMED[®]-classic Sortiment

REF 083468

Product Description, Intended use, Application

The beige adhesive dressing is made of a textile rigid fabric which is provided with a well-adhering zinc oxide adhesive which has a high ageing resistance. RUDAMED[®]-classic is used to cover smaller primary healing wounds. The non-adherent wound pad absorbs blood and exudate. Three plaster strips of different widths (4 cm, 6 cm, 8 cm) are packed in one plastic bag.

Composition

Cotton, polyethylene, viscose, zinc oxide adhesive

Contra indications

The product should not be used in the case of a known allergy against the material.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Legal Requirements

Medical product according to the directive MDD 93/42/EEC and the regulation MDR (EU 2017/745).

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary:	PE pouch
Secondary:	carton made of cellulose

Symbols used in labelling

Explanations at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and ISO 20417.

Storage and Transport

Dry and dustfree

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.