



RUDAMED[®]-Junior

REF 084106

Product Description, Intended use, Application

The carrier material of the adhesive dressing (6 cm x 10 cm) is made of a water-resistant, flexible synthetic printed with color designs for children. RUDAMED[®]-junior is used to cover smaller, primary healing wounds. The non-adherent wound pad absorbs blood and exsudate. Owing to its polyacrylate adhesive coating the product is suitable for sensitive skin.

Composition

Polyethylene, polypropylene, viscose and polyacrylate adhesive

Contra indications

The product should not be used in the case of a known allergy against the material.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Legal Requirements

Medical product according to the directive MDD 93/42/EEC and the regulation MDR (EU 2017/745).

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary:	PE bag
Secondary:	folding box made of cellulose
Tertiary:	carton made of cellulose

Symbols used in labelling

Explanations at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and ISO 20417.

Storage and Transport

Dry and dustfree

Single use device

Reusing a single use medical device can lead to microbiological danger. Re-processing, for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.