

NOBAMED Paul Danz AG

NOBADENT®

REF 149410

Product Description, Intended use, Application

The disposable product NOBADENT® is a white, soft, chlorine-free bleached cotton wool roll, size 1 cm x 4 cm. NOBADENT® is used for dental treatment in the dental field and serves to absorb saliva, secretion, and blood. The dental roll has a high absorption capacity, at the same time has a padding function and holding off function.

Composition

Cotton

Notes

The non-sterile cotton wool roll must be sterilized in accordance with a validated sterilization method before use on open wounds. The product can be sterilized with moist heat at 121°C/ 134°C, 2 to 3 bar, according to DIN EN 17665, radiation according to DIN EN ISO 11137 or ethylene oxide according to DIN EN ISO 11135.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulative Requirements, Common Standards

Medical device according to the MDD 93/42/EEC and MDR (EU) 2017/745.

Cotton wool according to BP or USP absorbent cotton

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging: folding box of

cellulose

Secondary packaging: carton made of

cellulose

Symbols used in labelling

Explanations at www.nobamed.com









Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and ISO 20417.

Storage and Transport

To be stored in a dry and dust-free environment

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.

Product Data Sheet: July 22, 2021 [REV 5] replaced: April 01, 2020 [Rev 4]