



NOBASURGMENT® ARTERIENKLEMME

REF 17-251

Product Description, Intended use, Application

Sterile single use artery clamp, Rochester, straight, for use in health care. The violet safety marking ensures clear visual differentiation from reusable instruments.

Composition

Stainless Steel

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulative Requirements, Common Standards

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging: paper-film packaging

Secondary packaging: carton made of cellulose

Symbols used in labelling

Explanations at www.nobamed.com



40-70% RH

Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417.

Storage and Transport

Dry and dustfree, protect from sunlight.

Sterile device

Before using a sterile product, visually inspect the packaging to ensure that it is intact.

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.