



**NOBAMED Paul Danz AG**

## NOBATISSUE<sup>®</sup>-weiß

REF 206546 (ab LOT 139152)

### Product Description, Intended use, Application

Non-sterile, **highly bleached** dressing cellulose, several individual and creped layers placed on top of each other, **size 37 cm x 57 cm**. They are absorbent, white, soft, and odourless. The highly bleached NOBATISSUE<sup>®</sup>-weiß products are mainly used in the health and care sector as absorbent material and for protection or padding purposes. NOBATISSUE<sup>®</sup>-weiß is for single use only.

### Composition

Highly bleached (oxygen bleach) dressing cellulose

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Regulative Requirements, Common Standards

Medical Device according to *MDR (EU) 2017/745*.

The cellulose layers comply with the requirements for dressing cellulose of the German Pharmacopoeia

The product does not contain dangerous toxic substances according to REACH.

### Packaging

Primary packaging: Carton made of cellulose

### Symbols used in labelling

Explanation at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.