

NOBAMED Paul Danz AG

# **NOBASLIP®-basic**

# **REF 310100**

# Product Description, Intended use, Application

The white, elastic net pants, **size M**, with **blue** colour coding at the waist to distinguish the sizes, is soft, well-fitting and skin-friendly. NOBASLIP<sup>®</sup>-basic serves for the fixation of medical devices such as absorbent pads or wound pads. The net pants are for single use.

## Composition

Polyester, elastane

## **Contra-Indications**

The product should not be used in the case of a known allergy against the material.

## **Notes**

Net pants which are too small will cause constrictions, if too big the fixation will be insufficient.

## **Incident reporting**

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

# Normative and Regulative

## **Requirements, Common Standards**

Medical device according to the MDD 93/42/EEC and MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

## Packaging

Primary packaging:

PE bag

Secondary packaging:

carton made of cellulose

## Symbols used in labelling Explanations at <u>www.nobamed.com</u>



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and ISO 20417.

### **Storage and Transport**

To be stored in a dry and dust-free environment

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.