

### NOBAMED Paul Danz AG

### **APPLIKATOR**

### **REF 320003**

# Product Description, Intended use, Application

The medical APPLICATOR is a cylindrical device (3 cm, for finger/ toe), consisting of two rings connected by longitudinal rods. The APPLICATOR is used for quick and easy application and placement of tubular bandages (such as NOBATRICOT®) over limbs during wound care.

### Composition

Stainless steel

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

# Normative and Regulative Requirements, Common Standards

Accessory for medical device tubular bandage according to the MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

#### **Packaging**

Primary packaging: PE-bag

Secondary packaging: carton made of

cellulose

### Symbols used in labelling

Explanations at www.nobamed.com







Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and DIN EN ISO 20417.

### **Storage and Transport**

To be stored in a dry and dust-free environment

#### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.