

## **NOBAMED Paul Danz AG**

## **NOBATRICOT®**

#### **REF 320080**

# **Product Description, Intended use, Application**

The white elastic tubular stockinet is 8 cm wide and 20 m long. It is made of low particle, hydrophilic polyester fibres and is hypoallergenic. NOBATRICOT is used to protect the skin as an underlayment for rigid and semi-rigid cast bandages. It can also be used as underlayment or cover for zinc paste bandages.

The tubular stockinet is also suited for the fixation of wound pads.

### Composition

Polyester

### **Contraindications**

The product should not be used in the case of a known allergy against the material.

#### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

# Normative and Regulative Requirements, Common Standards

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

### **Packaging**

Primary packaging: folding box

made of cellulose

Secondary packaging: carton made of cellu-

lose

### Symbols used in labelling

Explanation at www.nobamed.com







Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and ISO 20417.

### **Storage and Transport**

Dry and dustfree.

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

### **Disposal**

According to locally applicable legal regulations and standards of infection prophylaxis.

Product Data Sheet: August 20, 2021 [Rev6], replaces: April 01, 2020 [Rev5]