



## NOBATRICOT<sup>®</sup>

REF 320160

### Product Description, Intended use, Application

The white elastic tubular stockinet is 16 cm wide and 20 m long. It is made of low particle, hydrophilic polyester fibres and is hypoallergenic. NOBATRICOT is used to protect the skin as an underlayment for rigid and semi-rigid cast bandages. It can also be used as underlayment or cover for zinc paste bandages.

The tubular stockinet is also suited for the fixation of wound pads.

### Composition

Polyester

### Contraindications

The product should not be used in the case of a known allergy against the material.

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Regulative Requirements, Common Standards

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

### Packaging

Primary packaging: folding box  
made of cellulose

Secondary packaging: carton made of cellulose

### Symbols used in labelling

Explanation at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and ISO 20417.

### Storage and Transport

Dry and dustfree.

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.