

## **NOBAMED Paul Danz AG**

### **NOBATRICOT®**

### **REF 320240**

# Product Description, Intended use, Application

The white elastic tubular stockinet is 24 cm wide and 20 m long. It is made of low particle, hydrophilic polyester fibres and is hypoallergenic. NOBATRICOT is used to protect the skin as an underlayment for rigid and semirigid cast bandages. It can also be used as underlayment or cover for zinc paste bandages.

The tubular stockinet is also suited for the fixation of wound pads.

### Composition

Polyester

#### **Contraindications**

The product should not be used in the case of a known allergy against the material.

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

# Normative and Regulative Requirements, Common Standards

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

### **Packaging**

Primary packaging: folding box

made of cellulose

Secondary packaging: carton made of cellu-

lose

### Symbols used in labelling

Explanation at <u>www.nobamed.com</u>



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and ISO 20417.

### **Storage and Transport**

Dry and dustfree.

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

#### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.

Product Data Sheet: August 20, 2021 [Rev6], replaces: April 01, 2020 [Rev5]