



NOBALINGE®

REF 320260

Product Description, Intended use, Application

Ready-to-use elastic bandage for the foot. It is used for the secure and quick fixation of wound pads in cases of injuries. The novel conception of the bandage ensures a fast and safe application.

Composition

Polyamide

Contraindications

The product should not be used in the case of a known allergy against the material.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Legal Requirements

Medical product according to the directive MDD 93/42/EEC and the regulation MDR (EU 2017/745).

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary:	folding box made of cellulose
Secondary:	carton made of cellulose

Symbols used in labelling

Explanations at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1

Storage and Transport

Dry and dustfree

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.