



NOBAMED Paul Danz AG

NOBAFAST®

REF 330030

Product Description, Intended use, Application

The white, tight, and elastic tubular bandage, **size 1 (colour coded in red)**, is soft, conformable, and skin-friendly. The colour marking of the individual widths helps the user to quickly choose the right width. NOBAFAST®, being a well-fitting bandage, serves as a protection for the skin and is used underneath rigid and half rigid cast bandages. Due to the elastane integrated in the fabric, the bandage has a high elasticity and can therefore be applied wrinkle-free when used to fix wound pads. The rolled bandage is 10 m long. It is individually packed and for single use.

Composition

Viscose, elastane, polyamide, polyester

Contra-Indications

The product should not be used in the case of a known allergy against the material.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulative Requirements, Common Standards

Medical device according to the MDD 93/42/EEC and MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging: folding box
made of
cellulose

Secondary packaging: carton made of
cellulose

Symbols used in labelling

Explanations at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and ISO 20417.

Storage and Transport

To be stored in a dry and dust-free environment

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.