



## NOBAFAST®

REF 330050

### Product Description, Intended use, Application

The white, tight, and elastic tubular bandage, **size 2 (colour coded in green)**, is soft, conformable, and skin-friendly. The colour marking of the individual widths helps the user to quickly choose the right width. NOBAFAST®, being a well-fitting bandage, serves as a protection for the skin and is used underneath rigid and half rigid cast bandages. Due to the elastane integrated in the fabric, the bandage has a high elasticity and can therefore be applied wrinkle-free when used to fix wound pads. The rolled bandage is 10 m long. It is individually packed and for single use.

### Composition

Viscose, elastane, polyamide, polyester

### Contra-Indications

The product should not be used in the case of a known allergy against the material.

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Regulative Requirements, Common Standards

Medical device according to the MDD 93/42/EEC and MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

### Packaging

Primary packaging: folding box made of cellulose

Secondary packaging: carton made of cellulose

### Symbols used in labelling

Explanations at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and ISO 20417.

### Storage and Transport

To be stored in a dry and dust-free environment

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.