



## NOBAFROTT®

REF 350056

### Product Description, Intended use, Application

The disposable product NOBAFROTT® is an unbleached, elastic terry cloth tube with a **size of 6 cm x 5 m**. Synthetic threads are incorporated in the back area to achieve optimal elasticity. Due to the natural fibers, it is suitable for patients prone to allergies. NOBAFROTT® is used as an underlay to protect the skin and as padding for plaster and synthetic bandages. It can be individually cut to size, is easy to apply, soft, vapor- and air-permeable, adapts to the shape of the body and is therefore very comfortable to wear. The plaster or synthetic bandage can be applied directly over the terry cloth tube without additional padding material.

### Composition

Cotton, elastane, polyamide

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Regulative Requirements, Common Standards

Medical device according to the MDD 93/42/EEC and MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

### Packaging

Primary packaging:                    folding    box  
    made      of  
    cellulose

Secondary packaging:                carton made of  
    cellulose

### Symbols used in labelling

Explanations at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and DIN EN 1041.

### Storage and Transport

To be stored in a dry and dust-free environment

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.