



## NOBAFROTT®-Fertigverband

REF 350635

### Product Description, Intended use, Application

The disposable product NOBAFROTT®-Fertigverband, is an unbleached, elastic ready-to-use tubular bandage made of terry cloth for the **forearm**. Synthetic threads are incorporated in the back area to achieve optimum elasticity. Due to the natural fibers, it is suitable for patients who are prone to allergies. NOBAFROTT®-Fertigverband is used as an underlay for skin protection and as padding for plaster and plastic bandages. It is easy to apply, soft, supple, permeable to water vapor and air, with thumb inlay and thus offers a high degree of wearing comfort.

### Composition

Cotton, elastane, polyamide

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Regulative Requirements, Common Standards

Medical device according to the MDD 93/42/EEC and MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

### Packaging

Primary packaging:                    folding     box  
    made        of  
    cellulose

Secondary packaging:                carton made of  
    cellulose

### Symbols used in labelling

Explanations at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and DIN EN 1041.

### Storage and Transport

To be stored in a dry and dust-free environment

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.