

### **NOBAMED Paul Danz AG**

### **HALSKRAWATTE**

REF: 373106

# **Product Description, Intended use, Application**

The cervical collar (size: children) has an anatomical fit and a convenient hook and loop fastener. The cover offers a good wearing comfort. The cervical collar is used to immobilize the cervical spine.

The bandage is, fasteners closed, washable at 30° C in the case of domestic use. Due to its material composition, the product may not be treated chemothermally and may therefore not be used on different patients.

#### Composition

Cotton, polyurethan, polyamide

### **Contraindications**

The product should not be used in the case of a known allergy against the material

#### **Incident reporting**

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## Normative and Legal Requirements

Medical product according to directive MDD 93/42/EEC and the regulation MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

**Packaging** 

Primary: folding box

made of cellulose

Secondary: carton made of

cellulose

### Symbols used in labelling

Explanations at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1-and ISO 20417.

### **Storage and Transport**

Dry and dust free

### Single use device

Intended for use on one patient. The reprocessing of single use devices for reuse can decrease the product's performance significantly. Reusing a single use device can lead to microbiological danger.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.