



## HALSKRAWATTE

REF: 373108

### Product Description, Intended use, Application

The cervical collar (8 cm, small neck circumference, 33 cm - 38 cm) has an anatomical fit and a convenient hook and loop fastener. The cover offers a good wearing comfort. The cervical collar is used to immobilize the cervical spine.

The bandage is, fasteners closed, washable at 30° C in the case of domestic use. Due to its material composition, the product may not be treated chemothermally and may therefore not be used on different patients.

### Composition

Cotton, polyurethan, polyamide

### Contraindications

The product should not be used in the case of a known allergy against the material.

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Legal Requirements

Medical product according to directive MDD 93/42/EEC and the regulation MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

### Packaging

Primary: folding box  
made of  
cellulose

Secondary: carton made of  
cellulose

### Symbols used in labelling

Explanations at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and ISO 20417.

### Storage and Transport

Dry and dust free

### Single use device

Intended for use on one patient. The re-processing of single use devices for reuse can decrease the product's performance significantly. Reusing a single use device can lead to microbiological danger.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.