



## NOBADUR®-Verband

REF 400209

### Product Description, Intended Use, Application

NOBADUR®-Verband contains, according to the indication, two brown ultra-short stretch bandages NOBADUR, 8 or 10 cm wide and each 5 m long (stretched), with a stretchability of approx. 38%. If the bandage shall reach to the crotch, it can be completed with the appropriate NOBADUR® bandages. NOBADUR®-Verband is used for the application of a strongly compressing Pütter bandage for the compression of the extremities in phlebology or lymphology and to support and relief in the case of distortions and contusions in traumatology and sports medicine.

### Composition

Cotton, nylon

### Notes

The product has been subjected to a validated microbial reduction process and is therefore pre-sterilized.

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established

### Normative and Regulative Requirements, Common Standards

Medical device according to the MDD 93/42/EEC and MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

### Packaging

Primary packaging: folding box made of cellulose

Secondary packaging: carton made of cellulose

### Symbols used in labelling

Explanations at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and ISO 20417.

### Storage and Transport

To be stored in a dry and dust-free environment

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.