

# NOBAMED Paul Danz AG

# **NOBAFIX®**

### **REF 410425**

# **Product Description, Intended use, Application**

This fixation bandage with a width of 12 cm and a stretched length of 4 m is a rolled, lengthwise-elastic dressing retention bandage. Due to its weave, the bandage has a fine, smooth and tight structure. The elastic bandage is used to fix wound pads, especially over joints and body curves. As a damp bandage, NOBAFIX® is also used to shape and/or fix plaster splints and cast dressings. NOBAFIX® is for single use.

### Composition

Polyester

#### **Contraindications**

The product should not be used in the case of a known allergy against the material.

## **Incident reporting**

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

# Normative and Regulative Requirements, Common Standards

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

### **Packaging**

Primary packaging: folding box

made of cellulose

Secondary packaging: carton

made of cellulose

## Symbols used in labelling

Explanations at www.nobamed.com





Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN 1041.

### **Storage and Transport**

Dry and dustfree, protect from sunlight.

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.

Product Data Sheet: 03.05.2021 [REV 7], replaces: 11.05.2020 [REV 6]