



**NOBAMED Paul Danz AG**

# NOBATEX®

## Individually packed in foil

**REF 413400**

### Product Description, Intended use, Application

This white, woven bandage with a **width of 10cm** and a **stretched length of 4m** is a lengthwise-elastic dressing retention bandage. Due to its weave, the bandage has a curled or creped structure. The elastic bandage is used to fix wound pads in which its padding effect is advantageous and is particularly suitable over joints and body curves. As a damp bandage, NOBATEX® is also used to shape and/or fix plaster splints. NOBATEX® is for single use.

### Composition

Polyester, polyamide

### Contraindications

The product should not be used in the case of a known allergy against the material.

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Regulative Requirements, Common Standards

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

NOBATEX® is manufactured according to DIN 61634.

### Packaging

Primary packaging: PP foil

Secondary packaging: folding box made of cellulose

Tertiary packaging: carton made of cellulose

### Symbols used in labelling

Explanations at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN 1041.

### Storage and Transport

Dry and dustfree.

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.