



**NOBAMED Paul Danz AG**

# NOBAHAFT<sup>®</sup>-crepp

**REF 465204**

## **Product Description, Intended use, Application**

The cohesive elastic bandage, 4 cm x 20 m, consists of a fabric with crepped surface and does not contain latex. It adheres only to itself, but not to skin or clothing and is used for the fixation of wound pads.

## **Composition**

Polyamide, Viscose, latex-free polymer

## **Contraindications**

The product should not be used in the case of a known allergy against the material.

## **Note**

During application, cohesive bandages should only be stretched slightly as otherwise the arterial blood supply might be disturbed.

## **Incident reporting**

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## **Normative and Regulatory Requirements, Common Standards**

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

## **Packaging**

Primary packaging: folding box made of cellulose

Secondary packaging: Carton made of cellulose

## **Symbols used in labelling**

Explanations at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN 1041.

## **Storage and Transport**

Dry and dustfree.

## **Single use device**

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

## **Disposal**

According to locally applicable legal regulations and standards of infection prophylaxis.