



## NOBAHAFT<sup>®</sup>-crepp

REF 465212

### Product Description, Intended use, Application

The cohesive elastic bandage, 12 cm x 20 m, consists of a fabric with crepped surface and does not contain latex. It adheres only to itself, but not to skin or clothing and is used for the fixation of wound pads.

### Composition

Polyamide, Viscose, latex-free polymer

### Contraindications

The product should not be used in the case of a known allergy against the material.

### Note

During application, cohesive bandages should only be stretched slightly as otherwise the arterial blood supply might be disturbed.

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Regulative Requirements, Common Standards

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

### Packaging

Primary packaging: folding box made of cellulose

Secondary packaging: Carton made of cellulose

### Symbols used in labelling

Explanations at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN 1041.

### Storage and Transport

Dry and dustfree.

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.