

# **NOBAMED Paul Danz AG**

# **NOBA®-PRETAPE**

#### RFF 506007

# Product Description, Intended use, Application

NOBA®-Pretape (size: 7 cm x 27.5 m) is a white, open-pore and thin foam bandage. It is used to fix and compress padding material required for the application of rigid support bandages during preparation. It reduces the penetration of moisture. The bandage can also be used as an underlay under tape bandages to protect the skin. NOBA®-Pretape is easy to apply as it doesn't roll up during application.

## Composition

Polyurethane foam

#### Contraindication

The product should not be used in the case of a known allergy against the material.

#### **Incident reporting**

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## Normative and Regulative Requirements, Common Standards Medical Device according to MDD 93/42/EEC,

MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

# **Packaging**

Primary: folding box

made of cellulose

Secondary: carton made of

cellulose

## Symbols used in labelling

Explanation at www.nobamed.com





Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417

#### **Storage and Transport**

Dry and dustfree, protect from sunlight.

## Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

#### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.