

NOBAMED Paul Danz AG

# **RUHRSTERN®-HAFTPUFFER**

### **REF 570222**

## Product Description, Intended use, Application

RUHRSTERN<sup>®</sup>-HAFTPUFFER (crutch tip) is a three-part capsule for walking aids. It consists of a low-abrasion, elastic base and a rigid top with an included metal nozzle. This structure absorbs shocks and thus helps to protect the user's joints. The feature of a joint-function-like capsule supports the user's stable gait and safe standing, even if the walking aid is slanted when it hits the ground. The capsule is to be mounted onto a walking aid with an outside diameter of **22 mm**.

The capsule is a disposable product and is packed in pairs.

#### Composition

Polyurethane, polyvinyl chloride, metal

#### **Incident reporting**

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

#### **Normative and Regulative**

**Requirements, Common Standards** Accessory for a medical device according to *MDD* 93/42/EEC, *MDR* (EU) 2017/745.

RUHRSTERN<sup>®</sup>-HAFTPUFFER complies with the requirements of DIN EN ISO 24415-1 and of ISO/FDIS 24415-2 relating to durability, friction and stability.

The product does not contain dangerous toxic substances according to REACH.

#### Packaging

Primary packaging:	folding made cellulose	box of
Secondary packaging:	carton made of cellulose	

#### Symbols used in labelling

Explanation at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417

#### **Storage and Transport**

No special requirements

#### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

#### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.