

**NOBAMED** Paul Danz AG

# **NOBAPAD®-steril**

## **REF 601815**

## Product Description, Intended use, Application

Especially soft and stretchable sterile padding bandage which is made of binder-free synthetic wadding. The bandage can be used as preparatory padding for treatments with plaster, polyester or fiberglass bandages. This bandage is likewise suited for the padding of Cramer splints, Volkmann splints or other unpadded splints, or else as padding underneath compression bandages.

#### Composition

Polyester

#### Contraindications

The product should not be used in the case of a known allergy against the material.

#### **Incident reporting**

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

#### Normative and Regulative Requirements, Common Standards

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

Sterilization complies with DIN EN 11137.

Pa	cka	gi	ng	

Primary packaging:

Secondary packaging:

film-film packaging

> folding box made of cellulose

Tertiary packaging:

Carton made of cellulose

## Symbols used in labelling

Explanations at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and ISO 20417.

### **Storage and Transport**

Dry and dust free.

#### **Sterile device**

Before using a sterile product, visually inspect the packaging to ensure that it is intact.

#### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

#### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.