



## NOBA<sup>®</sup> transfer sheet

REF 620150

### Product Description, Intended use, Application

Disposable transfer sheet (75 cm x 190 cm) with four handles incorporated in the material on each side. It serves both as an absorbent underlay for stretchers, examination tables and beds and for lateral transfer of patients. The sheet creates a mechanical barrier between the patient and the surface of the bed etc., thus counteracting the spread of germs. The absorption volume according to ISO 11948-1 is appr. 3000 mL.

### Composition

Polypropylene, polyethylene, cellulose, superabsorber

### Contraindication

The product should not be used in the case of a known allergy against the material.

### Note

For patient up to 150 kg. For a weight of more than 150 kg up to a maximum of 230 kg, two sheets must be used on top of each other.

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Regulative Requirements, Common Standards

Medical Device according to MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

### Packaging

Primary:	PE bag
Secondary:	carton made of cellulose

### Symbols used in labelling

Explanations at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and DIN EN ISO 20417

### Storage and Transport

Dry and dustfree

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.