



NOBAMED Paul Danz AG

NOBA[®]-NASENSCHLITZTUCH (NOBA[®] nose slit cloth)

REF 651100

Product description, intended use, application

Face cloth for treatment tables in physiotherapy, dimensions: 40 cm x 30 cm, with Y-slot. The product is used to counteract the spread of germs and to protect the treatment table from sweat and cosmetic residues.

Composition

Airlaid cellulose

Contraindications

The product should not be used in the case of a known allergy against the material.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulative Requirements, Common Standards

Medical Device according to MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging: plastic bag

Secondary packaging: carton made of cellulose

Symbols used in labelling

Explanation at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and DIN EN ISO 20417

Storage and Transport

Dry and dustfree.

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.