

NOBAMED Paul Danz AG

# PFLEGEKITTEL

### **REF 664105**

## Product Description, Intendend use, Application

The protection gown, **size M**, is made of nonwoven. Its front and arms are impermeable to liquid, blood, and viruses. The yellow protection gown is open in the back and has knitted, elastic cuffs as well as ribbons for tiefastening in the waist and at the collar. The cuffs allow for a tight fit of disposable gloves. Due to the breathable material of the backside, the product is comfortable to wear. The yellow colour of the protection gown can be used for the assignment to specific sectors for example when nursing MRSA and VRSA patients. The use of the gown reduces the spread of germs and contaminations during nursing or medical activities. The gown is for single use.

#### Composition

Polypropylene (SPP), Polyethylene (front and arms)

Cuffs: mixture of polyester and cotton

#### **Contraindications**

The product should not be used in the case of a known allergy against the material.

#### **Incident reporting**

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

#### Normative and Regulative Requirements, Common Standards

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745. The product does not contain dangerous toxic substances according to REACH. The impermeability of the protection gown's front and arms against liquids, blood, and viruses has been proven in accordance with ASTM F 1670 "Synthetic Blood Test", ASTM F 1671 "Viral Penetration Test", ASTM F 739 "Test for Permeation of Liquids", ASTM F 903 "Test of resistance to penetration by liquids", AAMI PB 70- AATC 127 "Liquid barrier performance".

#### Packaging

Primary packaging: PE- bag

Secondary packaging: carton made of cellulose

#### Symbols used in labelling

Explanations at www.nobamed.com



**Storage and Transport** Dry and dustfree.

#### Single use device

Reusing aa single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

#### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.